


## Phospho-Cdc6-S54 Rabbit Mab For Wb, Elisa - Q99741

Item Number: CM0001942



### Introduction

Phospho-Cdc6-S54 Rabbit mAb (Q99741) validated for WB, ELISA. Detects human and mouse CDC6 phosphorylated at Ser54. Ideal for DNA replication and cell cycle checkpoint research. CamelBio one-stop IVD raw materials.

[Learn More](#)

| Parameter        | Value                                                                 |
|------------------|-----------------------------------------------------------------------|
| Product Name     | Phospho-Cdc6-S54 Rabbit mAb                                           |
| Remarks/Alias    | CDC18L; HsCDC6; MGORS5; HsCDC18; Phospho-Cdc6-S54                     |
| Species          | Human                                                                 |
| GeneID (Human)   | 990                                                                   |
| GeneID           | 990                                                                   |
| Immunogen        | A synthetic phosphorylated peptide around S54 of human Cdc6 (Q99741). |
| Source           | Rabbit                                                                |
| Category         | Phosphorylated Antibodies                                             |
| Application      | WB, ELISA                                                             |
| Cross Reactivity | Mouse                                                                 |
| SWISS            | Q99741                                                                |
| Protein Weight   | 63kDa                                                                 |
| Shipping         | Ice bag                                                               |