


Phospho-Cdk1-T14 Rabbit Mab For Wb, Elisa - P06493

Item Number: CM0001634



Introduction

Phospho-specific rabbit monoclonal antibody targeting CDK1 phosphorylated at Thr14. Validated for WB and ELISA in human and mouse samples. Detects a ~34 kDa band. Ideal for cell cycle, mitosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK1-T14 Rabbit mAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; Phospho-CDK1-T14
Species	Human
GeneID (Human)	983
GeneID	983
Immunogen	A synthetic phosphorylated peptide around T14 of human CDK1 (NP_001777.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06493
Protein Weight	34 kDa
Shipping	Ice bag