


Phospho-Cdk1-T161 Rabbit Pab - P06493

Item Number: CM0001656



Introduction

Phospho-CDK1-T161 Rabbit pAb specifically detects CDK1 phosphorylated at Thr161 in human, mouse, and rat samples. Validated for WB, IP, IF-P, IHC-P, and ELISA. Ideal for cell cycle and mitosis studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK1-T161 Rabbit pAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; Phospho-CDK1-T161
Species	Human
Gene ID	983
Immunogen	A synthetic phosphorylated peptide around T161 of human CDK1 (NP_001777.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag