


Phospho-Cdk1-T161+Cdk2/Cdk3-T160 Rabbit Monoclonal Antibody For Wb, Ip, If - P06493

Item Number: CM0000667



Introduction

High-specificity rabbit monoclonal antibody recognizing phospho-T161 of CDK1 and corresponding phospho-T160 of CDK2/CDK3. Validated for Western blot, immunoprecipitation, and immunofluorescence in human, mouse, and rat samples. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK1-T161+CDK2/CDK3-T160 Rabbit mAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; CDK1
Species	Human
GeneID (Human)	983
GeneID	983
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T161 of human CDK1 (NP_001777.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF-P, IHC-P, DB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag