


Phospho-Cdk1-Y15 Rabbit Mab For Wb, Elisa - P06493

Item Number: CM0001140



Introduction

Phospho-CDK1-Y15 Rabbit mAb for WB, ELISA. Detects phosphorylated CDK1 at Tyr15 in human samples. Ideal for cell cycle, mitosis, and DNA damage research.

[Learn More](#)

Specification	Detail
Product Name	Phospho-CDK1-Y15 Rabbit mAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; Phospho-CDK1-Y15
Species	Human
GeneID (Human)	983
GeneID	983
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y15 of human CDK1(NP_001777.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag