

Phospho-Cdk1/2-T14 Rabbit Monoclonal Antibody For Wb, Elisa - P06493 / P24941

Item Number: CM0000875



Introduction

Phospho-Cdk1/2-T14 Rabbit Monoclonal Antibody validated for WB and ELISA; detects phosphorylated T14 of human CDK1 and CDK2, cross-reacts with mouse and rat; suitable for cell cycle and mitosis studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Cdk1/2-T14 Rabbit mAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; Phospho-Cdk1/2-T14
Species	Human
Gene ID (Human)	983 / 1017
Gene ID	983 / 1017
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T14 of human CDK1 (P06493) / CDK2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493 / P24941
Protein Weight	34kDa
Shipping	Ice bag