


Phospho-Cdk2-Y15 Rabbit Mab For Wb, Elisa - P24941

Item Number: CM0000883



Introduction

Phospho-specific rabbit monoclonal antibody targeting human CDK2 phosphorylated at Tyr15. Validated for WB and ELISA, reactive with human and mouse. Ideal for cell cycle and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK2-Y15 Rabbit mAb
Remarks/Alias	CDKN2; p33(CDK2); Phospho-CDK2-Y15
Species	Human
GeneID (human)	1017
GeneID	1017
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around Y15 of human CDK2 (NP_001789.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24941
Protein Weight	34 kDa
Shipping	Ice bag