


Phospho-Cdk6-Y13 Rabbit Pab - Q00534

Item Number: CM0001833



Introduction

Phospho-CDK6-Y13 Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA. Detects phosphorylated Y13 of human CDK6; cross-reacts with mouse. Suitable for cell cycle and differentiation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK6-Y13 Rabbit pAb
Remarks/Alias	MCPH12; PLSTIRE; Phospho-CDK6-Y13
Species	Human
Gene ID	1021
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y13 of human CDK6 (NP_001250.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00534
Protein Weight	37kDa
Shipping	Ice bag