


Phospho-Cdk6-Y24 Rabbit Pab For Wb/If/Icc - Q00534

Item Number: CM0000413



Introduction

Phospho-CDK6-Y24 Rabbit pAb detects human CDK6 phosphorylated at Tyr24. Suitable for WB, IF/ICC. Validated in human, mouse. Ideal for cell cycle, differentiation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK6-Y24 Rabbit pAb
Remarks/Alias	MCPH12; PLSTIRE; Phospho-CDK6-Y24
Species	Human
GeneID (Human)	1021
GeneID	1021
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding Y24 of human CDK6
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human, Mouse
SWISS	Q00534
Protein Weight	37kDa
Shipping	Ice bag