


Phospho-Cdkn1B/P27Kip1-T187 Rabbit Pab For Wb - P46527

Item Number: CM0001583



Introduction

Phospho-specific rabbit polyclonal antibody targeting CDKN1B/p27KIP1 phosphorylated at T187; validated for Western blot in human samples. Ideal for cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDKN1B/p27KIP1-T187 Rabbit pAb
Remarks/Alias	KIP1; MEN4; CDKN4; MEN1B; P27KIP1; Phospho-CDKN1B/p27KIP1-T187
Species	Human
GeneID (human)	1027
GeneID	1027
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding T187 of human CDKN1B/p27KIP1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P46527
Protein Weight	22kDa
Shipping	Ice bag