

Phospho-Cebpb (T235) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P17676

Item Number: CM0000775



Introduction

Phospho-CEBPB (T235) rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat CEBPB. Ideal for transcription regulation, immune response, and adipogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CEBPB-T235 Rabbit mAb
Remarks/Alias	TCF5; IL6DBP; NF-IL6; C/EBP-beta
Species	Human
Gene ID (Human)	1051
Gene ID	1051
Immunogen	A synthetic phosphorylated peptide around T235 of human CEBPB (NP_005185.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17676
Protein Weight	36kDa
Shipping	Ice bag