

Phospho-Cfl1-S3 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P23528

Item Number: CM0000502



Introduction

Phospho-CFL1-S3 Rabbit pAb detects endogenous phosphorylated cofilin-1 at Ser3 in human, mouse, and rat. Suitable for WB, IHC-P, IF/ICC, ELISA. A reliable IVD tool for studying actin dynamics and cytoskeletal regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CFL1-S3 Rabbit pAb
Remarks/Alias	CFL; cofilin; HEL-S-15; Phospho-CFL1-S3
Species	Human
GeneID (Human)	1072
GeneID	1072
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S3 of human CFL1 (NP_005498.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23528
Protein Weight	19kDa
Shipping	Ice bag