

Phospho-Cgas (S420) Rabbit Polyclonal Antibody For Wb And Elisa - Q8C6L5

Item Number: CM0000665



Introduction

Phospho-CGAS-S420 Rabbit pAb detects CGAS phosphorylated at serine 420 in mouse, human, and rat samples. Validated for WB and ELISA. Ideal for innate immunity, cGAS-STING pathway, and DNA sensing research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-CGAS-S420 Rabbit pAb
Remarks/Alias	Mb21d1; E330016A19Rik; Phospho-CGAS-S420
Species	Mouse
GeneID (Human)	115004
GeneID	214763
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S420 of mouse CGAS (NP_775562.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8C6L5
Protein Weight	58kDa
Shipping	Ice bag