


Phospho-Chk1-S280 Rabbit Pab - O14757

Item Number: CM0000542



Introduction

Phospho-specific rabbit polyclonal antibody targeting Ser280 of human Chk1 (CHEK1) for WB, IHC-P, ELISA. Validated in human, mouse, rat. Ideal for DNA damage and cell cycle checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Chk1-S280 Rabbit pAb
Remarks/Alias	CHK1; Phospho-Chk1-S280
Species	Human
Gene ID (Human)	1111
Gene ID	1111
Immunogen	Synthetic peptide: a synthetic phosphorylated peptide around S280 of human Chk1 (NP_001265.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag