

Phospho-Chk2-S516 Rabbit Pab - O96017

Item Number: CM0001269



Introduction

Rabbit polyclonal antibody detecting endogenous levels of Chk2 only when phosphorylated at Ser516. Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for DNA damage response and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Chk2-S516 Rabbit pAb
Remarks/Alias	CDS1; CHK2; LFS2; RAD53; hCds1; HuCds1; PP1425; Phospho-Chk2-S516
Species	Human
GeneID	11200
Immunogen	A synthetic phosphorylated peptide around S516 of human Chk2 (NP_009125.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O96017
Protein Weight	61kDa
Shipping	Ice bag