


Phospho-Connexin 43-S368 Rabbit Polyclonal Antibody For IHC-P, Elisa - P17302

Item Number: CM0001217



Introduction

Phospho-specific rabbit pAb targeting human connexin-43 phosphorylated at Ser368. Validated for IHC-P and ELISA; cross-reacts with mouse, rat. Ideal for gap junction, cardiac, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Connexin 43-S368 Rabbit pAb
Remarks/Alias	HSS; CMDR; CX43; EKVP; GJAL; ODDD; AVSD3; EKVP3; HLHS1; PPKCA; Phospho-Connexin 43-S368
Species	Human
GeneID (Human)	2697
GeneID	2697
Immunogen	A synthetic phosphorylated peptide around S368 of human GJA1 (NP_000156.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P17302
Protein Weight	43kDa
Shipping	Ice bag