

## Phospho-Creb1-S129 Rabbit Pab - P16220-1

Item Number: CM0001507

### Introduction

Phospho-CREB1-S129 Rabbit pAb specifically detects CREB1 phosphorylated at Ser129 in human and rat samples. Validated for WB and ELISA, this antibody supports signal transduction and transcription regulation studies.

[Learn More](#)

| Parameter        | Value                                                         |
|------------------|---------------------------------------------------------------|
| Product Name     | Phospho-CREB1-S129 Rabbit pAb                                 |
| Remarks/Alias    | CREB; CREB-1; Phospho-CREB1-S129                              |
| Species          | Human                                                         |
| Gene ID (Human)  | 1385                                                          |
| Immunogen        | A synthetic phosphorylated peptide around S129 of human CREB1 |
| Source           | Rabbit                                                        |
| Category         | Phosphorylated Antibodies                                     |
| Application      | WB, ELISA                                                     |
| Cross Reactivity | Human, Rat                                                    |
| SWISS            | P16220-1                                                      |
| Protein Weight   | 36KDa                                                         |
| Shipping         | Ice bag                                                       |