

Phospho-Creb1-S133 Rabbit Monoclonal Antibody For Wb, Ihc-P, Fc (Intra), Elisa - P16220

Item Number: CM0000888



Introduction

Phospho-specific rabbit monoclonal antibody targeting CREB1 phosphorylated at Ser133. Validated applications include WB, IHC-P, intracellular flow cytometry, and ELISA with human, mouse, rat cross-reactivity. Suitable for cAMP signaling and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CREB1-S133 Rabbit mAb
Remarks/Alias	CREB; CREB-1
Species	Human
GeneID (human)	1385
GeneID	1385
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S133 of human CREB1 (NP_604391.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16220
Protein Weight	35kDa
Shipping	Ice bag