

## Phospho-Creb1-S133 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P16220

Item Number: CM0001280



### Introduction

Phospho-CREB1-S133 Rabbit pAb detects phosphorylated CREB1 at Ser133 in human, mouse, and rat. Suitable for WB, IHC-P, IF/ICC, and ELISA, enabling transcription and cAMP signaling research.

[Learn More](#)

| Parameter        | Value                                                                                          |
|------------------|------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-CREB1-S133 Rabbit pAb                                                                  |
| Remarks/Alias    | CREB; CREB-1; Phospho-CREB1-S133                                                               |
| Species          | Human                                                                                          |
| Gene ID (Human)  | 1385                                                                                           |
| Gene ID          | 1385                                                                                           |
| Immunogen        | Synthetic peptide; A synthetic phosphorylated peptide around S133 of human CREB1 (NP_604391.1) |
| Source           | Rabbit                                                                                         |
| Category         | Phosphorylated Antibodies                                                                      |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                              |
| SWISS            | P16220                                                                                         |
| Protein Weight   | 35kDa                                                                                          |
| Shipping         | Ice bag                                                                                        |