

Phospho-Creb1-S142 Rabbit Pab For Wb - P16220

Item Number: CM0000674



Introduction

Phospho-CREB1-S142 Rabbit pAb specifically detects CREB1 phosphorylated at Ser142 by Western blotting. Validated in human, mouse, and rat samples, it supports research on transcription regulation, circadian biology, and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CREB1-S142 Rabbit pAb
Remarks/Alias	CREB; CREB-1; Phospho-CREB1-S142
Species	Human
Gene ID (Human)	1385
Gene ID	1385
Immunogen	Synthetic peptide - A phospho specific peptide corresponding to residues surrounding S142 of human CREB1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P16220
Protein Weight	35kDa
Shipping	Ice bag