


Phospho-Crk P38-Y221 Rabbit Monoclonal Antibody For Wb, Elisa - P46108

Item Number: CM0001059



Introduction

Phospho-specific rabbit monoclonal antibody against Crk p38 at Y221. Validated for WB and ELISA, human cross-reactivity. Suitable for studying cell adhesion, migration, and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Crk p38-Y221 Rabbit mAb
Remarks/Alias	p38; CRKII
Species	Human
GeneID (Human)	1398
GeneID	1398
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y221 of human Crk p38 (NP_058431.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P46108
Protein Weight	34 kDa
Shipping	Ice bag