


Phospho-Cyclin D1 (T286) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24385

Item Number: CM0001049



Introduction

Rabbit polyclonal antibody against human Cyclin D1 phosphorylated at T286, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying G1/S transition and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Cyclin D1-T286 Rabbit pAb
Remarks/Alias	BCL1; PRAD1; U21B31; D11S287E; Phospho-Cyclin D1-T286
Species	Human
GeneID (Human)	595
GeneID	595
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T286 of human Cyclin D1 (NP_444284.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24385
Protein Weight	34kDa
Shipping	Ice bag