

Phospho-Cyclin E1-T395 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P24864

Item Number: CM0000776



Introduction

Phospho-specific rabbit polyclonal antibody targeting Cyclin E1 phosphorylated at T395. Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Cyclin E1-T395 Rabbit pAb
Remarks/Alias	CCNE; pCCNE1; Phospho-Cyclin E1-T395
Species	Human
GeneID (Human)	898
GeneID	898
Immunogen	A synthetic phosphorylated peptide around T395 of human Cyclin E1 (NP_001229.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24864
Protein Weight	47kDa
Shipping	Ice bag