


Phospho-Dna-Pkcs-S2056 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P78527

Item Number: CM0001805



Introduction

Rabbit polyclonal antibody specifically recognizing human Phospho-DNA-PKcs (S2056). Validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for DNA damage and repair signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-DNA-PKcs-S2056 Rabbit pAb
Remarks/Alias	HYRC; p350; DNAPK; DNPK1; HYRC1; IMD26; XRCC7; DNAPKc; DNA-PKC; DNA-PKcs; Phospho-DNA-PKcs-S2056
Species	Human
GeneID	5591
Immunogen	Synthetic phosphorylated peptide around S2056 of human DNA-PKcs (NP_008835.5)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P78527
Protein Weight	469kDa
Shipping	Ice bag