

Phospho-Dok2-Y299 Rabbit Pab - O60496

Item Number: CM0001560



Introduction

Phospho-DOK2-Y299 Rabbit pAb specifically detects human DOK2 phosphorylated at tyrosine 299. Validated for Western blot and immunofluorescence applications, this antibody is essential for studying cytokine and growth factor signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-DOK2-Y299 Rabbit pAb
Remarks/Alias	p56DOK; p56dok-2; Phospho-DOK2-Y299
Species	Human
Gene ID (Human)	9046
Gene ID	9046
Immunogen	A phospho specific peptide corresponding to residues surrounding Y299 of human DOK2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human
SWISS	O60496
Protein Weight	45kDa
Shipping	Ice bag