


Phospho-Drp1-S616 Rabbit Mab - 000429

Item Number: CM0001350



Introduction

Phospho-DRP1-S616 Rabbit mAb for WB, IP, ELISA applications, detecting endogenous levels of DRP1 only when phosphorylated at S616 in human, mouse, and rat samples. Suitable for mitochondrial dynamics research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-DRP1-S616 Rabbit mAb
Remarks/Alias	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV
Species	Human
GeneID (Human)	10059
GeneID	10059
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S616 of human DRP1 (NP_036192.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00429
Protein Weight	82kDa
Shipping	Ice bag