

Phospho-E-Cadherin-S838/S840 Rabbit Pab - P12830

Item Number: CM0001081



Introduction

Rabbit polyclonal antibody specific for human E-cadherin phosphorylated at S838/S840.

Validated for WB and ELISA. Essential for cell adhesion, epithelial cell mobility, and invasive suppressor studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-E-Cadherin-S838/S840 Rabbit pAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324; Phospho-E-Cadherin-S838/S840
Species	Human
Gene ID (Human)	999
Gene ID	999
Immunogen	A synthetic phosphorylated peptide around S838 & S840 of human CDH1 (NP_004351.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P12830
Protein Weight	97kDa
Shipping	Ice bag