


Phospho-Egfr-Y1045 Rabbit Pab For Wb, Elisa - P00533

Item Number: CM0001631



Introduction

Phospho-specific EGFR-Y1045 rabbit polyclonal antibody validated for WB and ELISA. Detects human EGFR phosphorylated at Tyr1045, a key residue in signal transduction. Ideal for studying EGFR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1045 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1045
Species	Human
Gene ID (Human)	1956
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Y1045 of human EGFR (NP_005219.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag