


Phospho-Egfr-Y1068 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P00533

Item Number: CM0001266



Introduction

Phospho-specific rabbit polyclonal antibody targeting EGFR phosphorylated at Y1068. Suitable for WB, IP, ELISA, with cross-reactivity to human and mouse. Ideal for studying EGFR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1068 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1068
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	A synthetic phosphorylated peptide around Y1068 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag