


Phospho-Egfr-Y1086 Rabbit Pab For Wb, Elisa - P00533

Item Number: CM0001252



Introduction

Phospho-EGFR-Y1086 Rabbit pAb is a rabbit polyclonal antibody targeting human EGFR phosphorylated at tyrosine 1086. Validated for WB and ELISA, it detects the 134 kDa receptor tyrosine kinase, supporting research into EGFR-driven signaling pathways and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1086 Rabbit pAb
Remarks/Alias	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1086
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around Y1086 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag