


Phospho-Egfr (Y1148) Rabbit Polyclonal Antibody For Wb, Elisa - P00533

Item Number: CM0001955



Introduction

Phospho-EGFR-Y1148 rabbit polyclonal antibody validated for WB and ELISA. Detects human EGFR phosphorylated at tyrosine 1148, enabling study of EGFR signaling pathways and receptor activation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1148 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1148
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y1148 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag