


Phospho-Egfr-Y1197 Rabbit Pab For Wb - P00533

Item Number: CM0001979



Introduction

Phospho-EGFR-Y1197 Rabbit Polyclonal Antibody detects human EGFR phosphorylated at tyrosine 1197 in Western blot. Ideal for studying EGFR signaling, cancer research, and drug resistance mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1197 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1197
Species	Human
Gene ID (Human)	1956
Immunogen	A phospho specific peptide corresponding to residues surrounding Y1197 of human EGFR
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag