

Phospho-Egfr-Y845 Rabbit Pab For Wb And Elisa - P00533

Item Number: CM0000560



Introduction

Rabbit polyclonal antibody targeting phosphorylated tyrosine 845 of human EGFR, suitable for Western blot and ELISA. Ideal for studies on EGFR signaling, kinase activity, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y845 Rabbit pAb
Remarks/Alias	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y845
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	A synthetic phosphorylated peptide around Y845 of human EGFR
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag