


Phospho-Egfr-Y992 Rabbit Pab For Wb And Elisa - P00533

Item Number: CM0000589



Introduction

Phospho-specific rabbit polyclonal antibody targeting EGFR phosphorylated at Y992. Suitable for WB and ELISA with human and mouse reactivity. Ideal for EGFR signaling and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y992 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y992
Species	Human
GeneID (human)	1956
GeneID	1956
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y992 of human EGFR (NP_005219.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag