


Phospho-Eif2A-S49 Rabbit Polyclonal Antibody For Wb, Ihc-P - P05198

Item Number: CM0001580



Introduction

Phospho-eIF2 α -S49 rabbit polyclonal antibody for Western blot and immunohistochemistry. Detects phosphorylated Ser49 of eIF2 α (EIF2S1) in human, mouse, rat. Ideal for integrated stress response and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF2 α -S49 Rabbit pAb
Remarks/Alias	EIF2; EIF-2; EIF2A; EIF-2A; EIF-2alpha; Phospho-eIF2 α -S49
Species	Human
Gene ID (Human)	1965
Immunogen	Synthetic peptide, a phospho specific peptide corresponding to residues surrounding S49 of human EIF2S1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P
Cross Reactivity	Human, Mouse, Rat
SWISS	P05198
Protein Weight	36kDa
Shipping	Ice bag