

Phospho-Eif4E-S209 Rabbit Pab - P06730

Item Number: CM0001582



Introduction

Phospho-eIF4E-S209 Rabbit pAb for Western blotting and ELISA, detecting human and mouse eIF4E phosphorylated at serine 209. Ideal for translation initiation, mRNA export, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4E-S209 Rabbit pAb
Remarks/Alias	CBP; EIF4F; AUTS19; EIF4E1; eIF-4E; EIF4EL1; Phospho-eIF4E-S209
Species	Human
Gene ID (Human)	1977
Gene ID	1977
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S209 of human eIF4E (NP_001959.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06730
Protein Weight	25kDa
Shipping	Ice bag