


Phospho-Eif4Ebp1 (S65) Rabbit Monoclonal Antibody For Wb, Elisa - Q13541

Item Number: CM0001051



Introduction

High-quality phospho-specific rabbit monoclonal antibody targeting eIF4EBP1 phosphorylated at Ser65, validated for Western blot and ELISA. Cross-reacts with human and mouse species. Ideal for translation regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4EBP1-S65 Rabbit mAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; Phospho-eIF4EBP1-S65
Species	Human
Gene ID (Human)	1978
Gene ID	1978
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S65 of human Phospho-eIF4EBP1-S65 (NP_004086.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag