


Phospho-Eif4Ebp1 (T37/T46) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13541

Item Number: CM0001319



Introduction

Rabbit polyclonal antibody against phosphorylated eIF4EBP1 (T37/T46) for WB, IHC-P, IF/ICC, ELISA. Validated in human, mouse, rat. Regulates translation initiation via mTORC1; key for cancer and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4EBP1-T37/46 Rabbit pAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; Phospho-eIF4EBP1-T37/46
Species	Human
Gene ID (Human)	1978
Gene ID	1978
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T37 & T46 of human eIF4EBP1 (NP_004086.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag