


Phospho-Eif4Ebp1-T37/T46 Rabbit Mab For Wb And Elisa - Q13541

Item Number: CM0000658



Introduction

Phospho-specific rabbit monoclonal antibody targeting eIF4EBP1 phosphorylation at T37/T46. Validated for WB, ELISA in human and rat samples. Essential for studying mTORC1 signaling and translation regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4EBP1-T37/T46 Rabbit mAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; Phospho-eIF4EBP1-T37/T46
Species	Human, Mouse, Rat
Gene ID (human)	1978
Immunogen	A synthetic phosphorylated peptide around T37&T46 of human EIF4EBP1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag