


Phospho-Eif4Ebp1-T45 Rabbit Pab For Western Blot - Q13541

Item Number: CM0001458



Introduction

Phospho-specific eIF4EBP1 (T45) rabbit polyclonal antibody for Western blot. Detects human eIF4EBP1 phosphorylated at Thr45, a key translational repressor in mTORC1 signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4EBP1-T45 Rabbit pAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; Phospho-eIF4EBP1-T45
Species	Human
Gene ID (Human)	1978
Gene ID	1978
Immunogen	Synthetic phospho-specific peptide corresponding to residues surrounding Thr45 of human eIF4EBP1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag