


Phospho-Eif4G1-S1108 Rabbit Pab For Wb, Elisa - Q04637

Item Number: CM0001570



Introduction

Rabbit polyclonal antibody specific for phospho-EIF4G1 (Ser1108) validated for WB and ELISA. Detects endogenous EIF4G1 in human, mouse, and rat samples, supporting translation initiation and stress response research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EIF4G1-S1108 Rabbit pAb
Remarks/Alias	P220; EIF4F; EIF4G; EIF4GI; PARK18; EIF-4G1; Phospho-EIF4G1-S1108
Species	Human
Gene ID (Human)	1981
Gene ID	1981
Immunogen	Synthetic peptide / A synthetic phosphorylated peptide around S1108 of human EIF4G1 (NP_001278086.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04637
Protein Weight	175kDa
Shipping	Ice bag