

Phospho-Enos-S1177 Rabbit Pab For Wb, Elisa - P29474

Item Number: CM0000985



Introduction

Rabbit polyclonal antibody specific for human eNOS phosphorylated at Ser1177. Validated for Western blot and ELISA, ideal for studying endothelial NO signaling, angiogenesis, and vascular regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eNOS-S1177 Rabbit pAb
Remarks/Alias	eNOS; ECNOS; Phospho-eNOS-S1177
Species	Human
Gene ID	4846
Immunogen	Synthetic phosphorylated peptide around S1177 of human eNOS (NP_000594.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P29474
Protein Weight	133kDa
Shipping	Ice bag