

Phospho-Epha2-S897 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P29317

Item Number: CM0001293



Introduction

Phospho-EPHA2 (S897) rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Detects human and mouse EPHA2, a receptor tyrosine kinase involved in cell adhesion, migration, and virus entry. Suitable for IVD research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EPHA2-S897 Rabbit mAb
Remarks/Alias	ECK; CTPA; ARCC2; CTPP1; CTRCT6; Phospho-EPHA2-S897
Species	Human
GeneID	1969
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S897 of human EPHA2 (NP_004422.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P29317
Protein Weight	108kDa
Shipping	Ice bag