

Phospho-Epha2-S897 Rabbit Pab For Wb, Elisa - P29317

Item Number: CM0001185



Introduction

Phospho-specific rabbit polyclonal antibody detecting EPHA2 phosphorylated at serine 897, validated for WB and ELISA in human, mouse, and rat. Ideal for studying EPHA2 signaling in cancer, development, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EPHA2-S897 Rabbit pAb
Remarks/Alias	ECK; CTPA; ARCC2; CTPP1; CTRCT6; Phospho-EPHA2-S897
Species	Human
Gene ID (Human)	1969
Gene ID	1969
Immunogen	A synthetic phosphorylated peptide around S897 of human EPHA2 (NP_004422.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29317
Protein Weight	108kDa
Shipping	Ice bag