

Phospho-Epha2-Y594 Rabbit Pab For Wb And Elisa - P29317

Item Number: CM0001610



Introduction

Rabbit polyclonal antibody specific for phosphorylated Y594 of human EPHA2 (Ephrin type-A receptor 2). Validated for western blot and ELISA; cross-reacts with human and rat. Suitable for receptor tyrosine kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EPHA2-Y594 Rabbit pAb
Remarks/Alias	ECK; CTPA; ARCC2; CTPP1; CTRCT6; Phospho-EPHA2-Y594
Species	Human
Gene ID	1969
Immunogen	A synthetic phosphorylated peptide around Y594 of human EphA2 (NP_004422.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P29317
Protein Weight	108kDa
Shipping	Ice bag