


Phospho-ErbB4/Her4-Y1284 Rabbit Pab - Q15303

Item Number: CM0000596



Introduction

Phospho-specific rabbit polyclonal antibody against ErbB4/HER4 at tyrosine 1284, validated in WB and ELISA. Detects human, mouse, and rat ErbB4. Ideal for studying receptor tyrosine kinase signaling in development and disease.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ErbB4/HER4-Y1284 Rabbit pAb
Remarks/Alias	HER4; ALS19; p180erbB4; Phospho-ErbB4/HER4-Y1284
Species	Human
Gene ID	2066
Immunogen	A synthetic phosphorylated peptide around Y1284 of human ErbB4/HER4 (NP_005226.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15303
Protein Weight	147kDa
Shipping	Ice bag