


Phospho-Erk1-T202 + Erk2-T185 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28482 / P27361

Item Number: CM0001598



Introduction

Phospho-specific rabbit polyclonal antibody recognizing ERK1 (phospho T202) and ERK2 (phospho T185). Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Essential for MAPK/ERK signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ERK1-T202 + ERK2-T185 Rabbit pAb
Species	Human
Gene ID	5594 / 5595
Immunogen	Synthetic phosphorylated peptide around T202 of human ERK1/ERK2 (NP_002737.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482 / P27361
Protein Weight	38kDa/40kDa/43kDa
Shipping	Ice bag