

Phospho-Erk1 (T202/Y204) + Erk2 (T185/Y187) Rabbit Mab For Wb, Ihc-P, Elisa - P27361/P28482

Item Number: CM0001809



Introduction

Phospho-specific rabbit monoclonal antibody against ERK1 (T202/Y204) and ERK2 (T185/Y187). Validated for Western blot, IHC paraffin, and ELISA in human, mouse, and rat tissues. 42/44kDa.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ERK1-T202/Y204 + ERK2-T185/Y187 Rabbit mAb
Remarks/Alias	ERK-1; ERK1; ERT2; HS44KDAP; HUMKER1A; P44ERK1; P44MAPK; PRKM3; p44-ERK1; p44-MAPK; Phospho-ERK1-T202/Y204 + ERK2-T185/Y187
Species	Human
GeneID (Human)	5594/5595
GeneID	5595/5594
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T202 & Y204 of human ERK1 (P27361).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27361/P28482
Protein Weight	42kDa/44kDa
Shipping	Ice bag