


## Phospho-Erk3-S189 Rabbit Polyclonal Antibody For Wb, Elisa - Q16659

Item Number: CM0000551



### Introduction

Phospho-specific rabbit polyclonal antibody targeting ERK3 (MAPK6) at Ser189. Validated for Western blot and ELISA applications; cross-reacts with human, mouse, and rat samples. Ideal for MAP kinase signaling studies.

[Learn More](#)

| Parameter        | Value                                                                      |
|------------------|----------------------------------------------------------------------------|
| Product Name     | Phospho-ERK3-S189 Rabbit pAb                                               |
| Remarks/Alias    | ERK3; PRKM6; p97MAPK; HsT17250; Phospho-ERK3-S189                          |
| Species          | Human                                                                      |
| Gene ID (Human)  | 5597                                                                       |
| Gene ID          | 5597                                                                       |
| Immunogen        | A synthetic phosphorylated peptide around S189 of human ERK3 (NP_002739.1) |
| Source           | Rabbit                                                                     |
| Category         | Phosphorylated Antibodies                                                  |
| Application      | WB, ELISA                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                          |
| SWISS            | Q16659                                                                     |
| Protein Weight   | 83kDa                                                                      |
| Shipping         | Ice bag                                                                    |