


Phospho-Esrα-S167 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P03372

Item Number: CM0001921



Introduction

Phospho-specific rabbit polyclonal antibody targeting Estrogen Receptor α (ER α /ESR1) at Ser167 phosphorylation. Validated for WB, IF/ICC, ELISA; reacts with human. Suitable for studying estrogen signaling, cancer biology, and nuclear receptor pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ESR α -S167 Rabbit pAb
Remarks/Alias	ER; ESR; Era; ESRA; ESTRR; NR3A1; Phospho-ESR α -S167
Species	Human
GeneID (Human)	2099
GeneID	2099
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S167 of human Estrogen Receptor alpha (ESR1) (NP_000116.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P03372
Protein Weight	66kDa
Shipping	Ice bag